

For collection, transport, maintenance  
and long term freeze storage



**Advanced  
Flocked Swabs**



**Virus, Chlamydia, Mycoplasma and Ureaplasma**

# UTM

UNIVERSAL TRANSPORT MEDIUM

• **Room Temperature Stable** •

**One transport - multiple tests - no sample rejection**

- Culture - all viruses, Chlamydia, mycoplasma and ureaplasma
- Enzyme Immuno Assays (EIA)
- PCR
- Nucleic acid amplification assays including *Chlamydia trachomatis* and *Neisseria gonorrhoeae* (Ct/Ng).

**Product stable at room temperature for 12 months**

**Organism viability maintained at room temperature and on ice**

**Choice of tube size, medium fill volume and product format**

Variety of patient specimen collection kits available comprising tubes of media and swabs or bulk packs of media tubes only.

**Advanced nylon flocked swabs increase sample release**

UTM is supplied with unique flocked swabs which increase sample release and sensitivity in tissue culture or assays.

**Precision snap molded  
breakpoint on swab shafts**

Plastic swab applicators supplied with UTM kits feature exclusive molded breakpoints.



**Capture cap docks and locks-in  
the swab stick**

UTM tube caps incorporate a unique internal design which functions like a docking system automatically capturing the end of the swab stick and securing it firmly into the cap.



**Internal conical shaped tube**

Tubes stand upright on the bench without the need for tube racks.

**Glass beads in each tube** assist the release of sample from the swab and mixing prior to analysis when used with a vortex mixer.



**Copan Italia S.p.A.**

Via F. Perotti, 10 • 25125 Brescia, Italy

Tel: (+39) 030 268 7211 • Fax: (+39) 030 268 7250

E-mail: [customer-care@copanitalia.com](mailto:customer-care@copanitalia.com) web: [www.copaninnovation.com](http://www.copaninnovation.com)



### COPAN UTM ORDERING INFORMATION

<b>BULK TUBE (no swabs)</b>		
<b>Cat. No.</b>	<b>Description</b>	<b>Packaging</b>
330C	16x100 mm tube, <b>3 ml</b> UTM medium (bulk pack)	6 x 50 kits
350C	12x80 mm mini tube, <b>1 ml</b> UTM medium (bulk pack)	6 x 50 kits

<b>TUBE + FLOCKED SWABS</b>		
<b>Cat. No.</b>	<b>Description</b>	<b>Packaging</b>
346C	16x100 mm tube, <b>3 ml</b> UTM medium, 1 Regular FLOCKED Swab, Molded Breakpoint	6 x 50 kits
307C	16x100 mm tube, <b>3 ml</b> UTM medium, 1 Minitip FLOCKED Swab, Molded Breakpoint	6 x 50 kits
305C	16x100 mm tube, <b>3 ml</b> UTM medium, 1 Flexible Nasopharyngeal FLOCKED Swab, Molded Breakpoint	6 x 50 kits
359C	12x80 mm mini tube, <b>1 ml</b> UTM medium, 1 Regular FLOCKED Swab, Molded Breakpoint	6 x 50 kits
361C	12x80 mm mini tube, <b>1 ml</b> UTM medium, 1 Minitip FLOCKED Swab, Molded Breakpoint	6 x 50 kits
360C	12x80 mm mini tube, <b>1 ml</b> UTM medium, 1 Flexible Nasopharyngeal FLOCKED Swab, Molded Breakpoint	6 x 50 kits